



Figure 6.1: Growth Hormone Release from Prototype 3 PLGA/SIS Scaffold


A. Schematic: GH Release from PLGA Scaffold




B. GH Concentration at Scaffold-Tissue Interface



C. Spatial GH Concentration Profiles




D. Cumulative GH Release Profile




# Supplementary Figure S1: Polymer Degradation and Diffusivity Evolution

## Polymer Degradation Profile



## Time-Dependent Diffusivity



## MODEL 5: BIOACTIVE AGENT DELIVERY FROM DEGRADING SCAFFOLD

### Growth Hormone Release from PLGA/SIS Scaffold - Simulation Results

---

#### SCAFFOLD PARAMETERS

- Outer Radius: 1.25 cm
- Inner Radius: 1.05 cm
- Wall Thickness: 0.20 cm
- Graft Length: 15.0 cm
- Scaffold Volume: 21.677 cm<sup>3</sup>

#### PLGA PROPERTIES

- Lactide:Glycolide Ratio: 75:25
- Degradation Rate Constant: 0.05 day<sup>-1</sup>
- Initial Molecular Weight: 100.0 kDa

#### GROWTH HORMONE LOADING

- Initial GH Loading: 0.15 mg
  - Initial Concentration: 0.0069 mg/cm<sup>3</sup> (6.92 µg/cm<sup>3</sup>)
- 

#### THERAPEUTIC WINDOW ANALYSIS

- Target Range: 10.0 - 50.0 ng/mL
- Maximum Interface Concentration: 127.32 ng/mL
- Time Above Minimum (10.0 ng/mL): 37.2 days
- Time Within Therapeutic Window: 32.7 days

#### RELEASE KINETICS

- Total Release at 30 days: 97.4%
- Remaining in Scaffold: 3.86 µg
- Average Release Rate: 3.25%/day

#### KORSMEYER-PEPPAS MODEL FIT

- Release Rate Constant (k): 0.0476
  - Release Exponent (n): 0.853
  - Release Mechanism: Anomalous (non-Fickian) transport
- 

#### FEASIBILITY ASSESSMENT

- ✓ FEASIBLE: Therapeutic concentration maintained for ≥3 weeks
- Prototype 3 with 0.15 mg GH loading is suitable for advancement.

The tri-phasic release profile (burst → diffusion → erosion) is evident from the simulation results, confirming the expected behavior of PLGA-based delivery systems.

---

#### DESIGN RECOMMENDATIONS

The modeling results support the following design considerations:

1. PLGA 75:25 provides controlled degradation kinetics
2. SIS coating effectively prevents luminal GH loss
3. Scaffold geometry (tubular) enables radial diffusion to tissue
4. Time-dependent diffusivity captures degradation-enhanced release

#### Next Steps:

- Validate model predictions with in vitro release studies
  - Conduct biocompatibility and stability assessments
  - Evaluate tissue penetration depth in ex vivo models
  - Optimize loading based on target therapeutic duration
-